



PROFESSIONAL GRADE DISINFECTANT & ODOR CONTROL

Now Available  
Classic &  
New  
Fresh Herbal  
Scent

**FOR USE WITH  
AIRX FRESHENAIRES  
& SCENTINELS**

**EVEN WITH  
COMPETITIVE  
CABINETS**



**RX 19E MICROCELL  
ODOR COUNTERACTANTS**

RX 19E microcells are the nearest thing to a liquid in solid form. They hold practically 100% of their calculable volume of *AIRICIDE® Odor Counteractant*, yet they won't spill or drip like liquids. Nor will they skin or crust over as do gels. These space age microcells hold the *AIRICIDE®* until exposed to the air, then gradually, but steadily, release it by capillary action to change the size and shape of odorous molecules so it is

recognized as something different and more pleasant than it would have been.

Odors float in the air in the form of molecules that when sniffed in are identified by minute hair-like antenna...The microcell blocks are small and compact (2½x4x½) allowing the use of multiple blocks in a single circulator or wall cabinet and adapters are available allow use of 19E with just about any competitive cabinet.



Block is easily removed from poly bag without even touching block.



Use RX 19Es with Airx Freshenaire to control intense foul odors.



Use adapters in existing and competitive cabinets.



Scentinel Jr. controls release thru adjustable fins. No battery or fans.

**AIRX LABORATORIES, 1640 DELMAR DRIVE, FOLCROFT, PA. 19032-0037**  
PHONE (610) 534-8900 • FAX (610) 534-8912 • Email - sales@airxinfo.com • www.airxinfo.com

# AIRICIDE<sup>®</sup>, RX 19E MICROCELLS & CABINETS



## 1. AIRICIDE<sup>®</sup>

Not a mask, not a desensitizer, AIRICIDE<sup>®</sup> is formulated to eliminate foul odors by changing the geometrics of odorous molecules in the air, so when they reach cilia receptors in the nose they are not the bad odors they would have been.

## 2. MICROCELLS

A space age, rigid polymer foam that holds nominally 100% of its dimensional (2½" x 4" x ½") capacity. Then gradually and constantly releases its vapors into the air to neutralize foul odors from cancer, burns, human incontinence, animals, decay, bacterial growth, fire, smoke, liquor, beer, mildew, etc.

### SPECIFICATIONS

TYPE.....	compressed microcells
ODOR COUNTERACTANT.....	pure AIRICIDE <sup>®</sup>
RELEASE OF COUNTERACTANT.....	approx. 60-90 days
OLFACTORY DESENSITIZERS.....	none
PETROLEUM DISTILLATES.....	none
TOXIC.....	no
EFFECT ON METALS.....	non-corrosive
COLOR.....	white, non-staining
FREEZE-THAW STABILITY.....	no degradation
HEAT STABILITY.....	no degradation
SHELF LIFE.....	min. 1 year, closed
DRYING CHARACTERISTICS.....	no shrinking, no skinning
VERSATILITY.....	neutralizes odors from bacterial growth, oxidation reduction, chemicals, etc.

## 3. THE AIRX FRESHENAIRE

A portable case (13"x16"x4") with handle. Built-in extension cord and 110v electric fan to circulate AIRICIDE<sup>®</sup>, and cover to conserve microcells when not being used.

## 4. MICROCELL ADAPTERS

These foam pieces allow conversion of any deodorizer cabinet to the use of RX 19E Microcells. Competitive cabinets are more effective with AIRICIDE<sup>®</sup>, and adapters can hold two microcells. There's no need to discard one that still has life when installing a new block.

## 5. THE SCENTINEL WALL CABINET

Brushed metallic cabinet with computerized mechanics operated by a single 'D' battery. Light on front signals low battery. Can be programmed to turn off when room is not occupied. Use with 1,2 or 3 microcells, depending on need.

## 6. THE SCENTINEL JR.

The simple solution to many odor problems. Holds a single RX 19 microcell. No battery to replace or be stolen. No fan – Fins on side of cabinet allow adjustment for individual needs.

DISTRIBUTED BY: